### MONITORING

### **BLOCK AND PASTE**

MONITORING RODENT ACTIVITIES

Thanks to their attractive formulations, LIPHATECH monitoring baits (paste or block) are indispensable to track rodent activities, movements or presence in areas where conventional baits cannot be used.

Their enhanced formulations make them easy to change to baits when it is necessary for a quick consumption.

#### INTEGRATED BAITING STRATEGY AS PLACEBO BAITS

Associated with AEGIS Trap Rat and AEGIS Mouse Trap, placebo baits will show the same advantages than other LIPHATECH baits with high palatability that will guide rodents directly in the traps.

	PALLETIZATION (80x120)
Bucket of 4 kg of 20g block	80 buckets
Bucket of 5 kg of 10g paste	80 buckets



### **ACCESSORIES**

### **AEGIS STICK & GO TETHERING CABLE**

Bait stations securing cables – Easy to use and quick to install.

Strong, flexible 40cm and 1m long.


Reliably adheres to all clean surface – with the double sided tape.

	PALLETIZATION (80x120)
Carton of 15 sachets of 10 cables of 40cm	216 cartons
Carton of 15 sachets of 10 cables of 1m	216 cartons

#### **AEGIS HOOKS**

Solid hook used for attaching blocks on the rod - Ideal for sewer's use.

	PALLETIZATION (80x120)
Carton of 100 plastic hooks	360 cartons
Carton of 500 metal hooks	108 cartons



GENERATION DIFFERNIALONE

AEGIS®
SECURED BAIT STATIONS

MAKI

ROZOL° CHLOROPHACINONE

GENERAL CONDITIONS OF USE OF THIS DOCUMENT (hereinafter the "Document")

Edition February 2017. Replaces and supersedes any and all prior editions. LIPHATECH reserves the right to modify this Document at any time. The reader/user of this Document is supposed to verify, before any use, that the edition available to him is the last updated one. (contact@liphatech.fr).

This Document is not an official or contractual document. It shall not be considered as a commercial offer or a warranty of any level of efficacy of the products. The information contained herein shall not be considered as specifications or recommendations of use of the products. This information is based on the data existing at the date of the edition of this Document and only indicative. Any total or partial use, diffusion, or reproduction of this Document is subject to the prior written consent of LIPHATECH and must be done in

All logos, tradenames, slogans, trademarks mentioned in this Document are protected by intellectual property rights. Any right of use is subject to the prior written consent of their respective owners.

The sale of LIPHATECH's products and material mentioned in this document are governed by the General Conditions of Sale of LIPHATECH. Any purchaser is supposed to verify the G.C.S. in force at the date of the order by contacting LIPHATECH or visiting www.liphatech.fr website.

**LIPHATECH** - S.A.S au capital de 2.800.000€ - Bonnel - CS10005 - 47480 PONT DU CASSE - FRANCE - 442 688 206 R.C.S. AGEN Tél. : +33 05 53 69 35 70 - E-mail : contact@liphatech.fr

WWW.LIPHATECH.FR

# PROFESSIONAL AND CONSUMER RANGES RODENT CONTROL





- GROUPE DE SANGOSSE

Use biocides safely. Always read the label and product information before use.

AEGIS® is a Registered Trademark of LIPHATECH SAS

## **AEGIS BAIT STATIONS**

- // Aegis stations are closed and secured with a key.
- // Protect baits against dust and humidity.
- Protect non target animals, humans, kids and pets, from bait consumption.
- Long lasting of baits without damages for.
- Long lasting palatability for increased efficiency.
- Possibility to store and fix baits (feeder, rod, sticks) to prevent dispersibility.



Universal Key



### **AEGIS® RAT COMPACT**

Two positioning options for different baiting plans. Possible to store baits in dedicated areas. Horizontal or vertical bait pinning - Specific feeder for grains. Dimensions: 330 x 210 x 100 mm.

Suitable for majority of sites and control of both roof rat or brown rat.

	PALLETIZATION (100x120)
Carton of 6 bait stations	54 cartons

### **TUNNEL STATION**

Optimized size for cable pathways or beams baiting - Ideal for the control of roof rats. Two options for bait pinning - Dimensions: 330 x 110 x 90 mm.

Carton of 24 bait stations	20 cartons	
No.		
-		



### **AEGIS** MOUSE

Ideal size for mice control - Adapted to mice habits for multiple baiting points. Optimized size toughness ratio - Dimensions : 125 x 95 x 37 mm.

	PALLETIZATION (100x120)	
Carton of 20 bait stations	240 cartons	
Carton of 100 bait stations	48 cartons	3
	Part and City	

### **MECHANICAL SOLUTIONS**

// Quick results // Easy visual control // Alternative to bait aversion situations

### **AEGIS** RATTRAP

EASY AND SAFE RESET WITHOUT OPENING THE BOX With limited risks of injury - No neophobic effect if associated with AEGIS Compact.

Possible to attract with all kind of food source. Hygienic - Trap cannot be moved by rodents. Protected and temper resistant in every weather conditions. Associated with placebo bait as attractant.



ITES WHERE CHEMICAL SOLUTIONS RE NOT THE APPROPRIATE ANSWER

To protect a site containing a lot of alternative food sources. Situations where end users want a proof of efficiency. Site with low infestation - Maintained control all year round.









### AEGIS MOUSETRAP

AFE FOR PESTS AND CHILDREN

Easy and fun to use and handle - Reusable and discreet. Hygienic without contact with dead mice - Possibility to rod any solid food or placebo paste as attractant. Baiting from the underside with no pinching risk. Triggered status can be checked from a distance - Bodies concealed.



RE NOT THE APPROPRIATE ANSWER

To protect a site containing a lot of alternative food sources. Situations where end users require a proof of efficiency. Site with low infestation - Maintained control all year round. Ideal for both professionals and consumers.



	PALLETIZATION (80x120)
Carton of 20	150 cartons
Carton of 12 (individual pack only)	147 cartons

